Acid Detergent Lignin Procedure

<u>Materials:</u>

Polypropylene digestion tubes and glass rods Buchner funnels Filtering flasks Aluminum weighing pan Filtering flasks Glass filter papers (5.5 cm) Desiccator

Reagents:

Sulfuric Acid (72% by weight) Solution:

72% H₂SO₄ (add 2451 g sulfuric acid (96%) (VWR Cat # 2876-46, FW 98.08, CAS # 7664-93-9) to 817 g dd H₂O) standardized to specific gravity of 1.634 at 20 °C.

Procedure:

- 1. Place paper with residue from ADF procedure into a labeled 125 ml polypropylene digestion tube. Push to bottom with glass stirring rod. Save ADF pan.
- 2. Wearing safety glasses and gloves add 40 ml 72% H₂SO₄. Use glass rod to break up sample. Stir well. Leave glass rod in tube.
- 3. Stir once every hour. Digest 4 hours.
- 4. Filter contents onto glass filter papers under light vacuum.
- 5. Rinse tube and glass rod onto filter with hot dd water. Wash residue until free of acid (at least 250 ml).
- 6. Carefully transfer filter from funnel into aluminum weighing pan. Dry 8 hours or overnight in drying oven (100 °C).
- 7. Remove from oven and cool in desiccator. Record weight.
- 8. Ash dried residue in muffle furnace 8 hours or overnight (470 °C).
- 9. Remove slightly cooled samples and cool in desiccator. Record weight of ashed residue.

Calculations:

	А	В	С	D	E	F	G	Н
Sample	Sample Weight	DM	Before Ash	After Ash	% Lignin dmb	Lignin Mean	sd	% CV
#2	0.3387	$0.9700 \\ 0.9700$	1.4628 1.4369	1.4591 1.4334	1.126	1.15	0.033	2.89
	$E = \left[\frac{C-D}{A \times B}\right] \times 100$			1.4354	$1.126 = \left[\frac{1.4628 - 1.4591}{0.3387 \times 0.9700}\right] \times 100$			

To calculate F, G, and H see statistical analysis procedure.

Reference

Goering, H. K. and P. J. Van Soest. 1970. Forage fiber analyses. Agr. Handbook No 379, ARS, USDA.

Fiber (acid detergent) and Lignin (H₂SO₄) in Animal Feed, AOAC Official Method 973.18, Official Methods of Analysis, 18th Ed, Revision 3, 2010